

DEHS Disposal Guide

Glass Trash



Lab personnel must tape box closed (any tape is acceptable), mark top of the box "GLASS TRASH" and place near (not in) regular trash can.

Box cannot exceed 20 lbs. (9 kg)



If box is not pre-marked to convey hazard, must write or label box.

Pipettes & Pipette tips...



Because they readily puncture waste and trash bags and could be an injurious hazard potential to waste handlers, non-contaminated pipettes and pipette tips should be disposed in glass disposal boxes*, not in regular trash receptacles. All pipettes and pipette tips should not have any free liquids when discarded into a biohazard container or glass trash box.

** Alternate wording is acceptable as long as it describes the non-contaminated waste contents.*

