

Heater Probe in Upper Endoscopy

(Inject “pillow” before SB burn)

Lesion	Probe	Pressure	Energy	Pulses/ Site
PUD/ Dieulafoy	10 Fr	Very Firm	30 J	4
M-W Tear	7-10 Fr	Moderate	20 J	3
Angiodys plasia	7-10 Fr	Light	15 J	2

Heater Probe in Colonoscopy

(Inject “Pillow” before burn; tattoo after)

Lesion	Probe	Pressure	Energy	Pulses/ Site
Ulcer	10 Fr	Moderate	15 J	2
Stalk	10 Fr	Moderate	15-20 J	2
Diverticuli	10 Fr	Moderate	15 J	2
Cancer	10 Fr	Moderate	20 J	2
Angiodys plasia	7-10 Fr	Light	10 J	1

BICAP in Upper Endoscopy

(Inject “pillow” before SB burn)

Lesion	Probe	Pressure	Energy	Time/ Site
PUD/ Dieulafoy	10 Fr	Very Firm	20 W	10-14 sec
M-W Tear	7-10 Fr	Moderate	20 W	4 sec
Angiodys plasia	7-10 Fr	Light	15 W	2 sec

BICAP in Colonoscopy

(Inject “Pillow” before burn; tattoo after)

Lesion	Probe	Pressure	Energy	Time/ Site
Ulcer	10 Fr	Moderate	20 W	2 sec
Stalk	10 Fr	Moderate	20 W	2 sec
Diverticuli	10 Fr	Moderate	20 W	2 sec
Cancer	10 Fr	Moderate	20 W	2 sec
Angiodysplasia	7-10 Fr	Light	15 W	1 sec